

Anti-Ftcd Monoclonal Antibody For Wb, Elisa - 095954

Item Number: CM0007488

Introduction

High-specificity recombinant rabbit monoclonal antibody against human FTCD (LCHC1), validated for WB and ELISA in human, mouse, and rat. Ideal for folate metabolism and cytoskeletal studies.

[Learn More](#)

Parameter	Value
Product Name	FTCD Rabbit PolymAb®
Remarks/Alias	LCHC1
Species	Human
GenelD (Human)	10841
GenelD	10841
Immunogen	Recombinant protein of human FTCD
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95954
Protein Weight	59kDa
Shipping	Ice bag