

Anti-Ftcd Monoclonal Antibody For Wb, Elisa - 095954

Item Number: CM0006650

Introduction

FTCD Rabbit Monoclonal Antibody validated for Western blot and ELISA. Recognizes human, mouse, and rat FTCD with a predicted molecular weight of 59 kDa. Ideal for histidine metabolism and folate pathway research.

[Learn More](#)

Parameter	Value
Product Name	FTCD Rabbit mAb
Remarks/Alias	LCHC1; FTCD
Species	Human
GenelD (human)	10841
GenelD	10841
Immunogen	Recombinant protein of human FTCD
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95954
Protein Weight	59 kDa
Shipping	Ice bag