

Anti-Ftcd Polyclonal Antibody For Wb And Elisa - 095954

Item Number: CM0024893

Introduction

Rabbit polyclonal antibody targeting human FTCD (formiminotransferase-cyclodeaminase), raised against a C-terminal recombinant fragment. Validated for WB and ELISA, cross-reacts with mouse and rat. Research use only.

[Learn More](#)

Parameter	Value
Product Name	FTCD Rabbit pAb
Remarks/Alias	LCHC1; FTCD
Species	Human
GeneID (Human)	10841
GeneID	10841
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 412-541 of human FTCD (NP_006648.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095954
Protein Weight	59kDa
Shipping	Ice bag