

Anti-Fto Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9C0B1

Item Number: CM0002770



Introduction

KO Validated anti-FTO rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Specific to human, mouse, rat FTO (~58 kDa). Ideal for epitranscriptomics, obesity, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] FTO Rabbit mAb
Remarks/Alias	GDFD; ALKBH9; BMIQ14; TO
Species	Human
GeneID (human)	79068
GeneID	79068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 327-504 of human FTO (NP_001073901.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0B1
Protein Weight	58kDa
Shipping	Ice bag