

Anti-Fto Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9C0B1

Item Number: CM0009646



Introduction

High-performance FTO rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA. Reacts with human FTO, a key RNA demethylase involved in fat mass, obesity, and leukemia. Ideal for epigenetics and metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	FTO Rabbit mAb
Remarks/Alias	GDFD; ALKBH9; BMIQ14; FTO
Species	Human
GeneID (human)	79068
GeneID	79068
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FTO (Q9C0B1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q9C0B1
Protein Weight	58kDa
Shipping	Ice bag