

Anti-Fto Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9C0B1

Item Number: CM0018709



Introduction

FTO Rabbit Polyclonal Antibody detects human, mouse, and rat FTO (alpha-ketoglutarate-dependent dioxygenase). Validated for WB, IP, ELISA. Targets the fat mass and obesity-associated RNA demethylase for epitranscriptomics and adipogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	FTO Rabbit pAb
Remarks/Alias	GDFD; ALKBH9; BMIQ14; FTO
Species	Human
GeneID (Human)	79068
GeneID	79068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 327-504 of human FTO (NP_001073901.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C0B1
Protein Weight	58kDa
Shipping	Ice bag