


Anti-Ftsj3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8ly81

Item Number: CM0029538



Introduction

Rabbit polyclonal antibody against human FTSJ3, a pre-rRNA 2'-O-ribose RNA methyltransferase. Validated for WB and ELISA. Ideal for research on ribosome biogenesis, rRNA processing, and HIV-1 innate immune evasion.

[Learn More](#)

Parameter	Value
Product Name	FTSJ3 Rabbit pAb
Remarks/Alias	SPB1; EPCS3; FTSJ3
Species	Human
GeneID (Human)	117246
GeneID	117246
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 608-847 of human FTSJ3 (NP_060117.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IY81
Protein Weight	97kDa
Shipping	Ice bag