


Anti-Fubp1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q96Ae4

Item Number: CM0025657



Introduction

Rabbit polyclonal antibody against human FUBP1 (FUSE-binding protein 1), validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse, rat. ~68kDa. For research use.

[Learn More](#)

Parameter	Value
Product Name	FUBP1 Rabbit pAb
Remarks/Alias	FBP; FUBP; hDH V; FUBP1
Species	Human
GeneID (Human)	8880
GeneID	8880
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-644 of human FUBP1 (NP_003893.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AE4
Protein Weight	68kDa
Shipping	Ice bag