

Anti-Fubp3 Polyclonal Antibody For Wb, Elisa - Q96I24

Item Number: CM0029252

Introduction

High-quality FUBP3 rabbit polyclonal antibody for WB and ELISA, detecting human and mouse far upstream element-binding protein 3. This nuclear DNA-binding protein may activate gene expression, making it a valuable tool for transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	FUBP3 Rabbit pAb
Remarks/Alias	FBP3; FUBP3
Species	Human
GeneID (Human)	8939
GeneID	8939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-260 of human FUBP3 (NP_003925.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96I24
Protein Weight	62kDa
Shipping	Ice bag