

Anti-Fukutin Monoclonal Antibody For Wb, Ihc-P, Elisa - 075072

Item Number: CM0003610

Introduction

Rabbit monoclonal antibody against human Fukutin (FKTN), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for dystroglycanopathy, congenital muscular dystrophy, and ribitol-phosphate transferase research.

[Learn More](#)

Parameter	Value
Product Name	Fukutin Rabbit mAb
Remarks/Alias	FCMD; CMD1X; LGMD2M; MDDGA4; MDDGB4; MDDGC4; LGMDR13; Fukutin
Species	Human
GeneID (Human)	2218
GeneID	2218
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 362-461 of human Fukutin (075072).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075072
Protein Weight	54kDa
Shipping	Ice bag