


Anti-Fundc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8lvp5

Item Number: CM0019396



Introduction

Rabbit polyclonal antibody against human FUNDC1, validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Targets a mitochondrial outer membrane protein involved in mitophagy and MAM formation. Ideal for autophagy and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	FUNDC1 Rabbit pAb
Remarks/Alias	FUNDC1
Species	Human
GeneID (human)	139341
GeneID	139341
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-155 of human FUNDC1 (NP_776155.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IVP5
Protein Weight	17kDa
Shipping	Ice bag