


Anti-Fundc1 Polyclonal Antibody For Wb - Q8lvp5

Item Number: CM0021612



Introduction

Rabbit polyclonal antibody against human FUNDC1, validated for Western blot (WB) and ELISA. Peptide immunogen from amino acids 50-155. Suitable for mitophagy, mitochondrial dynamics, and MAMs research.

[Learn More](#)

Parameter	Value
Product Name	FUNDC1 Rabbit pAb
Remarks/Alias	FUNDC1
Species	Human
GeneID (Human)	139341
GeneID	139341
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-155 of human FUNDC1 (NP_776155.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IVP5
Protein Weight	17kDa
Shipping	Ice bag