


## Anti-Fundc1 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - Q8lvp5

Item Number: CM0011310



### Introduction

Recombinant anti-FUNDC1 rabbit monoclonal antibody validated for WB, IP, IHC-P, and ELISA. Detects human, mouse, and rat FUNDC1 (17 kDa). Suitable for mitophagy, mitochondrial dynamics, and angiogenesis research.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | FUNDC1 Rabbit mAb                                                                                    |
| Remarks/Alias    | Fundc1                                                                                               |
| Species          | Human                                                                                                |
| GeneID           | 139341                                                                                               |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FUNDC1 (NP_776155.1) |
| Source           | Rabbit                                                                                               |
| Category         | Monoclonal Antibodies                                                                                |
| Application      | WB, IP, IHC-P, ELISA                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                    |
| SWISS            | Q8IVP5                                                                                               |
| Protein Weight   | 17kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |