

Anti-Furin Monoclonal Antibody For Wb And Elisa - P09958

Item Number: CM0010922



Introduction

Anti-Furin rabbit monoclonal antibody for WB and ELISA. Human and mouse reactive, targeting the proprotein convertase Furin. Suitable for research and diagnostic applications.

[Learn More](#)

Parameter	Specification
Product Name	Furin Rabbit mAb
Remarks/Alias	FUR; PACE; SPC1; PCSK3; Furin
Species	Human
Gene ID (Human)	5045
Gene ID	5045
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Furin (P09958).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09958
Protein Weight	87kDa
Shipping	Ice bag