

Anti-Furin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P09958

Item Number: CM0025667



Introduction

Rabbit polyclonal antibody against human Furin (PCSK3/PACE), suitable for Western blot, IF/ICC, and ELISA. Recognizes human, mouse, and rat Furin. Applications include research on TGF-beta signaling and host-pathogen interactions.

[Learn More](#)

Parameter	Value
Product Name	Furin Rabbit pAb
Remarks/Alias	FUR; PACE; SPC1; PCSK3; Furin
Species	Human
GeneID (Human)	5045
GeneID	5045
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Furin (NP_002560.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09958
Protein Weight	87kDa
Shipping	Ice bag