

Anti-Fus (Ko Validated) Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P35637

Item Number: CM0029112



Introduction

High-quality rabbit polyclonal antibody against human RNA-binding protein FUS (TLS), validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying RNA processing, DNA repair, and ALS-associated mechanisms.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] FUS Rabbit pAb
Remarks/Alias	TLS; ALS6; ETM4; FUS1; POMP75; altFUS; HNRNPP2; US
Species	Human
Gene ID	2521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 297-526 of human FUS (NP_004951.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35637
Protein Weight	53kDa
Shipping	Ice bag