

Anti-Fus Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P35637

Item Number: CM0018980



Introduction

FUS rabbit polyclonal antibody validated for WB, IHC-P, IP, ELISA; cross-reacts with human, mouse, rat. Supports research on RNA-binding protein FUS in neurodegenerative diseases and cancer.

[Learn More](#)

Parameter	Value
Product Name	FUS Rabbit pAb
Remarks/Alias	TLS; ALS6; ETM4; FUS1; POMP75; altFUS; HNRNPP2; FUS
Species	Human
Gene ID (Human)	2521
Gene ID	2521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 297-526 of human FUS (NP_004951.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35637
Protein Weight	53kDa
Shipping	Ice bag