

Anti-Fut11 Rabbit Polyclonal Antibody For Wb And Elisa - Q495W5

Item Number: CM0023629

Introduction

FUT11 rabbit polyclonal antibody, raised against recombinant human FUT11 protein, validated for WB and ELISA. Suitable for O-fucosylation research and glycobiology studies. Reacts with human FUT11 protein (~56 kDa).

[Learn More](#)

Parameter	Value
Product Name	FUT11 Rabbit pAb
Remarks/Alias	FUCTXI; FUT11
Species	Human
GenelD (Human)	170384
GenelD	170384
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q495W5
Protein Weight	56kDa
Shipping	Ice bag