

Anti-Fut2 Rabbit Polyclonal Antibody For Wb, Elisa - Q10981

Item Number: CM0030788

Introduction

Rabbit polyclonal antibody targeting human FUT2 (secretor factor), suitable for Western blot and ELISA. Validated for research use in glycosylation and blood group antigen studies.

[Learn More](#)

Parameter	Value
Product Name	FUT2 Rabbit pAb
Remarks/Alias	SE; Se2; sej; SEC2; B12QTL1; FUT2
Species	Human
GeneID (human)	2524
GeneID	2524
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q10981
Protein Weight	39kDa
Shipping	Ice bag