

Anti-Fut4 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P22083

Item Number: CM0024003

Introduction

Rabbit polyclonal antibody targeting human FUT4 (CD15/LeX). Suitable for WB, IF-P, IHC-P, and ELISA; cross-reacts with mouse and rat. Explore fucosyltransferase biology in leukocyte migration and development.

[Learn More](#)

Parameter	Specification
Product Name	FUT4 Rabbit pAb
Remarks/Alias	LeX, CD15, ELFT, FCT3A, FUTIV, SSEA-1, FUC-TIV, FUT4
Species	Human
GeneID (Human)	2526
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 186-331 of human FUT4 (NP_002024.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22083
Protein Weight	59kDa
Shipping	Ice bag