

Anti-Fut6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P51993

Item Number: CM0025284

Introduction

Rabbit polyclonal antibody targeting human fucosyltransferase 6 (FUT6/Fuc-TVI). Validated for WB, IHC-P, and ELISA. Designed for glycobiology research, including sialyl Lewis X biosynthesis and related biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	FUT6 Rabbit pAb
Remarks/Alias	FT1A; FCT3A; Fuc-TVI; FucT-VI; FUT6
Species	Human
Gene ID (Human)	2528
Gene ID	2528
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P51993
Protein Weight	42kDa
Shipping	Ice bag