

Anti-Fut7 Polyclonal Antibody For Wb, Elisa - Q11130

Item Number: CM0027514

Introduction

Rabbit anti-human FUT7 polyclonal antibody (FucT-VII) for WB, ELISA. Reacts with human, mouse, rat. Investigates fucosyltransferase 7 and sialyl Lewis X-mediated leukocyte adhesion. Q11130.

[Learn More](#)

Parameter	Value
Product Name	FUT7 Rabbit pAb
Remarks/Alias	FucT-VII; FUT7
Species	Human
GeneID (human)	2529
GeneID	2529
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human FUT7 (NP_004470.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q11130
Protein Weight	39kDa
Shipping	Ice bag