


Anti-Fxr1 Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - P51114

Item Number: CM0009064



Introduction

Rabbit monoclonal anti-FXR1 antibody for WB, IP, IF/ICC, ELISA. Detects endogenous FXR1 (~70 kDa) in human, mouse, rat. RNA-binding protein regulating mRNA translation/stability in neurogenesis, muscle, spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	FXR1 Rabbit mAb
Remarks/Alias	FXR1P; CMYP9A; CMYP9B; MYOPMIL; MYORIBF; FXR1
Species	Human
GeneID	8087
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human FXR1 (P51114).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51114
Protein Weight	70kDa
Shipping	Ice bag