

Anti-Fxr2 Rabbit Polyclonal Antibody For Wb, Elisa - P51116

Item Number: CM0017381



Introduction

Rabbit polyclonal antibody targeting FXR2, validated for Western blot (WB) and ELISA in human, mouse, and rat samples. Targets the RNA-binding protein FXR2, which regulates mRNA translation and is required for hippocampal neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	FXR2 Rabbit pAb
Remarks/Alias	FXR2P; FMR1L2; FXR2
Species	Human
Gene ID (Human)	9513
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 550-650 of human FXR2 (NP_004851.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51116
Protein Weight	74kDa
Shipping	Ice bag