


Anti-Fxr/Nr1H4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Ri1

Item Number: CM0006924



Introduction

High-quality rabbit monoclonal antibody targeting human FXR/NR1H4, validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for bile acid receptor signaling and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	FXR/NR1H4 Rabbit mAb
Remarks/Alias	BAR; FXR; HRR1; HRR-1; PFIC5; RIP14; FXR/NR1H4
Species	Human
Gene ID (Human)	9971
Gene ID	9971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human FXR/NR1H4 (NP_001193922.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Ri1
Protein Weight	48 kDa/54 kDa/55 kDa/56 kDa
Shipping	Ice bag