


Anti-Fxyd1 (Phospholemman) Rabbit Polyclonal Antibody For Wb, Elisa - 000168

Item Number: CM0019507



Introduction

Rabbit polyclonal antibody against human FXYD1/Phospholemman, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studies of ion transport regulation and NKA activity.

[Learn More](#)

Parameter	Value
Product Name	FXYD1 Rabbit pAb
Remarks/Alias	PLM; FXYD1
Species	Human
GeneID (Human)	5348
GeneID	5348
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-92 of human FXYD1 (NP_005022.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	000168
Protein Weight	10kDa
Shipping	Ice bag