

Anti-Fxyd3 Polyclonal Antibody For Ihc-P, Elisa - Q14802

Item Number: CM0019357

Introduction

High-quality FXYD3 polyclonal antibody validated for IHC-P and ELISA. Detects human FXYD3/Mat-8, a sodium/potassium-ATPase regulator, ideal for ion transport and cancer biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	FXYD3 Rabbit pAb
Remarks/Alias	MAT8; PLML; FXYD3
Species	Human
GeneID (Human)	5349
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-87 of human FXYD3 (NP_005962.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q14802
Protein Weight	9kDa
Shipping	Ice bag