

Anti-Fxyd6 Polyclonal Antibody For Wb, Elisa - Q9H0Q3

Item Number: CM0017723



Introduction

Rabbit polyclonal antibody against human FXYD6 (phosphohippolin), suitable for WB and ELISA. Cross-reacts with mouse and rat. Targets the sodium/potassium ATPase regulator for ion transport research.

[Learn More](#)

Parameter	Value
Product Name	FXYD6 Rabbit pAb
Remarks/Alias	FXYD6
Species	Human
Gene ID	53826
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human FXYD6 (NP_071286.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H0Q3
Protein Weight	11kDa
Shipping	Ice bag