

Anti-Fxyd7 Polyclonal Antibody For Wb, Elisa - P58549

Item Number: CM0016844

Introduction

Rabbit polyclonal antibody against human FXYD7, validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying sodium/potassium ATPase regulation and ion transport.

[Learn More](#)

Parameter	Value
Product Name	FXYD7 Rabbit pAb
Remarks/Alias	FXYD7
Species	Human
GeneID (human)	53822
GeneID	53822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human FXYD7 (NP_071289.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P58549
Protein Weight	9kDa
Shipping	Ice bag