


Anti-Fyb Polyclonal Antibody For Wb, Elisa - O15117

Item Number: CM0016252



Introduction

Rabbit polyclonal antibody against human FYB (FYN-binding protein 1), validated for WB and ELISA. Cross-reacts with human and mouse. Ideal for studying T-cell signaling and actin cytoskeleton remodeling.

[Learn More](#)

Parameter	Value
Product Name	FYB Rabbit pAb
Remarks/Alias	FYB; ADAP; THC3; PRO0823; SLAP130; SLAP-130
Species	Human
Gene ID (Human)	2533
Gene ID	2533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 678-829 of human FYB (NP_001456.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15117
Protein Weight	85kDa
Shipping	Ice bag