

Anti-Fyn Rabbit Monoclonal Antibody For Wb, Elisa - P06241

Item Number: CM0007901



Introduction

Highly validated rabbit monoclonal antibody against human Fyn (P06241), specific for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for tyrosine kinase signaling studies in immune response, axon guidance, and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	Fyn Rabbit mAb
Remarks/Alias	SLK, SYN, p59-FYN, Fyn
Species	Human
Gene ID (Human)	2534
Gene ID	2534
Immunogen	A synthetic peptide corresponding to amino acids 1-100 of human Fyn (P06241)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06241
Protein Weight	61kDa
Shipping	Ice bag