

## Anti-Fzd3 Polyclonal Antibody For Wb, Elisa - Q9Npg1

Item Number: CM0016105



### Introduction

Rabbit polyclonal antibody targeting human Frizzled-3 (FZD3), validated for WB and ELISA applications. Ideal for studying Wnt signaling, neurogenesis, and planar cell polarity. Cross-reacts with rat.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | FZD3 Rabbit pAb                                                                                                  |
| Remarks/Alias    | Fz-3; FZD3                                                                                                       |
| Species          | Human                                                                                                            |
| GeneID           | 7976                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 23-205 of human FZD3 (NP_059108.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Rat                                                                                                              |
| SWISS            | Q9NPG1                                                                                                           |
| Protein Weight   | 76kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |