

Anti-Fzd6 Polyclonal Antibody For Wb, Elisa - O60353

Item Number: CM0015146



Introduction

Rabbit polyclonal antibody targeting human Frizzled-6 (FZD6), suitable for WB and ELISA; cross-reacts with mouse and rat. Useful for Wnt signaling and planar cell polarity studies.

[Learn More](#)

Parameter	Specification
Product Name	FZD6 Rabbit pAb
Remarks/Alias	FZ6; FZ-6; HFZ6; NDNC1; NDNC10; FZD6
Species	Human
GeneID (Human)	8323
GeneID	8323
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-201 of human FZD6 (NP_003497.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60353
Protein Weight	79kDa
Shipping	Ice bag