

Anti-Fzd7 Rabbit Polyclonal Antibody For Wb, Elisa - 075084

Item Number: CM0030112



Introduction

High-quality rabbit polyclonal antibody against human Frizzled-7 (FZD7). Validated for WB and ELISA, detects human and mouse FZD7 (predicted 64 kDa). Ideal for Wnt signaling research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	FZD7 Rabbit pAb
Remarks/Alias	FzE3; FZD7
Species	Human
GeneID (Human)	8324
GeneID	8324
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human FZD7 (NP_003498.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075084
Protein Weight	64kDa
Shipping	Ice bag