

Anti-Fzd9 Rabbit Polyclonal Antibody For Wb And Elisa - O00144

Item Number: CM0030508

Introduction

Rabbit polyclonal antibody against human Frizzled-9 (FZD9, CD349), validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat samples. Ideal for Wnt signaling and neurodevelopment studies.

[Learn More](#)

Parameter	Value
Product Name	FZD9 Rabbit pAb
Remarks/Alias	FZD3, CD349, FZD9
Species	Human
Gene ID (Human)	8326
Gene ID	8326
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 469-591 of human FZD9 (NP_003499.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00144
Protein Weight	64kDa
Shipping	Ice bag