

Anti-G0S2 Polyclonal Antibody For Wb, Elisa - P27469

Item Number: CM0017287

Introduction

Rabbit polyclonal antibody targeting human G0S2, a mitochondrial protein that promotes apoptosis by binding BCL2. Validated for Western blot and ELISA, it is ideal for studies of cell death regulation and G0/G1 switch mechanisms.

[Learn More](#)

Parameter	Value
Product Name	G0S2 Rabbit pAb
Remarks/Alias	G0S2
Species	Human
Gene ID (Human)	50486
Gene ID	50486
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human G0S2 (NP_056529.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P27469
Protein Weight	11kDa
Shipping	Ice bag