

Anti-G0S2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P27469

Item Number: CM0029307

Introduction

Rabbit polyclonal antibody targeting human G0S2 protein, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. G0S2 promotes apoptosis by binding BCL2 and localizes to mitochondria. Ideal for apoptosis research and diagnostic assays.

[Learn More](#)

Parameter	Value
Product Name	G0S2 Rabbit pAb
Remarks/Alias	G0S2
Species	Human
GeneID (Human)	50486
GeneID	50486
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human G0S2 (NP_056529.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27469
Protein Weight	11kDa
Shipping	Ice bag