

## Anti-G2E3 Rabbit Polyclonal Antibody For Wb And Elisa - Q7L622

Item Number: CM0017357

### Introduction

Rabbit polyclonal antibody targeting human G2E3, an E3 ubiquitin ligase for histone H3K14 monoubiquitination. Validated for WB and ELISA; cross-reacts with human, mouse, rat. Ideal for chromatin and cell cycle studies.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | G2E3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | PHF7B; KIAA1333; G2E3                                                                                           |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 55632                                                                                                           |
| Gene ID          | 55632                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human G2E3 (NP_060239.2) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | Q7L622                                                                                                          |
| Protein Weight   | 81kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |