

Anti-G3Bp1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13283

Item Number: CM0023507



Introduction

Rabbit polyclonal antibody targeting human G3BP1, suitable for WB, IHC-P, IF/ICC, ELISA. Recognizes protein from human, mouse, rat. Essential tool for investigating stress granule dynamics, innate immunity, and nucleic acid helicase functions.

[Learn More](#)

Parameter	Value
Product Name	G3BP1 Rabbit pAb
Remarks/Alias	G3BP; HDH-VIII; G3BP1
Species	Human
Gene ID (Human)	10146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 131-330 of human G3BP1 (NP_938405.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13283
Protein Weight	52kDa
Shipping	Ice bag