

Anti-G3Bp1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13283

Item Number: CM0005001



Introduction

Rabbit monoclonal antibody targeting human G3BP1 (Ras GTPase-activating protein-binding protein 1). Validated for Western blotting, immunofluorescence/immunocytochemistry, immunoprecipitation, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for stress granule and innate immunity research.

[Learn More](#)

Parameter	Specification
Product Name	G3BP1 Rabbit mAb
Remarks/Alias	G3BP; HDH-VIII; G3BP1
Species	Human
GeneID (human)	10146
GeneID	10146
Immunogen	Synthetic Peptide corresponding to a sequence within amino acids 350-450 of human G3BP1 (Q13283)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13283
Protein Weight	52kDa
Shipping	Ice bag