

Anti-G3BP1 Rabbit Monoclonal Antibody [Ko Validated] For Wb, If/Icc, Ihc-P, Elisa - Q13283

Item Number: CM0004864

Introduction

High-quality rabbit monoclonal antibody against human G3BP1, validated for Western blot, IF/ICC, IHC-P, and ELISA. Recognizes human, mouse, and rat G3BP1. Suitable for stress granule research, innate immunity, and helicase studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] G3BP Rabbit mAb
Remarks/Alias	G3BP; HDH-VIII
Species	Human
GeneID (Human)	10146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-323 of human G3BP1 (NP_938405.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13283
Protein Weight	52 kDa/14 kDa
Shipping	Ice bag