

Anti-G3BP2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Un86

Item Number: CM0030626



Introduction

G3BP2 Rabbit pAb (polyclonal) validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets stress granule scaffold protein (54kDa) critical for antiviral signaling. Ideal for RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	G3BP2 Rabbit pAb
Remarks/Alias	G3BP2
Species	Human
Gene ID (Human)	9908
Gene ID	9908
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-280 of human G3BP2 (NP_987100.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UN86
Protein Weight	54kDa
Shipping	Ice bag