


Anti-G6Pc Rabbit Polyclonal Antibody For Wb, Elisa - P35575

Item Number: CM0019767



Introduction

Rabbit polyclonal antibody against G6PC (Glucose-6-phosphatase catalytic subunit 1). Validated for WB and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for studying glycogenolysis, gluconeogenesis, and glycogen storage disease.

[Learn More](#)

Parameter	Value
Product Name	G6PC Rabbit pAb
Remarks/Alias	G6PC; G6PT; GSD1; GSD1a; G6Pase
Species	Human
Gene ID (Human)	2538
Gene ID	2538
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35575
Protein Weight	40kDa
Shipping	Ice bag