

Anti-G6Pd Polyclonal Antibody For Wb, If/Icc, Elisa - P11413

Item Number: CM0017705



Introduction

High-quality rabbit polyclonal antibody targeting G6PD (G6PD1), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat species. Suitable for investigating glucose-6-phosphate dehydrogenase function in metabolic pathways, hemolytic anemia, and oxidative stress.

[Learn More](#)

Parameter	Value
Product Name	G6PD Rabbit pAb
Remarks/Alias	G6PD1; G6PD
Species	Human
GeneID	2539
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 195-515 of human G6PD (NP_001035810.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11413
Protein Weight	59kDa
Shipping	Ice bag