

Anti-G6Pd Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11413

Item Number: CM0011775



Introduction

High-affinity rabbit monoclonal antibody against human G6PD, validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat G6PD. Suitable for metabolic pathway and redox biology studies.

[Learn More](#)

Parameter	Value
Product Name	G6PD Rabbit mAb
Remarks/Alias	G6PD1
Species	Human
GeneID (Human)	2539
GeneID	2539
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 416-515 of human G6PD (NP_001035810.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11413
Protein Weight	60kDa
Shipping	Ice bag