


Anti-Gaa Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P10253

Item Number: CM0028260



Introduction

Rabbit polyclonal antibody targeting human lysosomal alpha-glucosidase (GAA), validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Essential for glycogen degradation; relevant to glycogen storage disease research.

[Learn More](#)

Parameter	Value
Product Name	GAA Rabbit pAb
Remarks/Alias	LYAG; GAA
Species	Human
GeneID (Human)	2548
GeneID	2548
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human GAA (NP_000143.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10253
Protein Weight	105kDa
Shipping	Ice bag