


Anti-Gab1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13480

Item Number: CM0026237



Introduction

High-quality rabbit polyclonal antibody against human GAB1, validated for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for studying GRB2-associated binding protein 1 in signaling cascades.

[Learn More](#)

Parameter	Value
Product Name	GAB1 Rabbit pAb
Remarks/Alias	DFNB26; GAB1
Species	Human
GeneID (Human)	2549
GeneID	2549
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 594-694 of human GAB1 (NP_002030.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13480
Protein Weight	77kDa
Shipping	Ice bag