

Anti-Gaba A Receptor Beta 1 (Gabrb1) Monoclonal Antibody For Wb, Elisa - P18505

Item Number: CM0003136



Introduction

Anti-GABA A Receptor beta 1 (GABRB1) Rabbit Monoclonal Antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for inhibitory neurotransmission and epilepsy research.

[Learn More](#)

Parameter	Value
Product Name	GABA A Receptor beta 1 (GABRB1) Rabbit mAb
Remarks/Alias	DEE45; EIEE45; GABA A Receptor beta 1 (GABRB1)
Species	Human
GeneID (Human)	2560
GeneID	2560
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of human GABA A Receptor beta 1 (GABRB1) (NP_000803.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18505
Protein Weight	54kDa
Shipping	Ice bag