


Anti-Gaba A Receptor Beta 2 (Gabrb2) Rabbit Polyclonal Antibody For Wb, Elisa - P47870

Item Number: CM0021808



Introduction

Rabbit polyclonal antibody specifically recognizing GABA A Receptor beta 2 (GABRB2). Suitable for WB and ELISA. Cross-reacts with human and mouse. A valuable tool for research on inhibitory neurotransmission and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	GABA A Receptor beta 2 (GABRB2) Rabbit pAb
Remarks/Alias	DEE92; ICEE2; GABA A Receptor beta 2 (GABRB2)
Species	Human
Gene ID (Human)	2561
Gene ID	2561
Immunogen	Recombinant fusion protein containing amino acids 26-244 of human GABRB2 (NP_000804.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47870
Protein Weight	59kDa
Shipping	Ice bag