


Anti-Gabarap Monoclonal Antibody For Wb, Elisa - O95166

Item Number: CM0009468



Introduction

Recombinant rabbit monoclonal antibody targeting human GABARAP, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for autophagy and intracellular transport research.

[Learn More](#)

Parameter	Value
Product Name	GABARAP Rabbit mAb
Remarks/Alias	MM46; ATG8A; GABARAP-a; GABARAP
Species	Human
GeneID (Human)	11337
GeneID	11337
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GABARAP (O95166)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95166
Protein Weight	14kDa
Shipping	Ice bag